

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. :RP10769

Lot No. :E1922326C

Time Processed :11:59:57

Month-Day-Year Processed :03/21/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
10.00	Pumps	B.Conc	35
35.00	Pumps	B.Conc	95
37.00	Pumps	B.Conc	95
37.01	Pumps	B.Conc	5
45.00	Pumps	B.Conc	5
45.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	21.175	73369	9617	0.504
2	24.350	410302	41010	2.820
3	24.717	52600	7103	0.361
4	25.092	13843002	1036523	95.127
5	25.742	172850	15526	1.188
Total		14552124	1109779	100.000