

Sample Name:KRAS WT Peptide (VVVGAGGVGK)
Cat No. : RP30235
Lot No. : E1975124A
Time Processed :15:34:56
Month-Day-Year Processed :01/15/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.658	21214	3062	0.180
2	6.750	14939	3063	0.126
3	6.975	11677363	799348	98.843
4	7.400	51793	9641	0.438
5	13.708	48784	4879	0.413
Total		11814093	819993	100.000