


Mass Spectrum



Sample Information		Equipment : ZJ21010035	
Month-Day Processed :	01/14/24	Interface Bias :	+4.5 kV
Time Processed :	11:00:24	Drying Gas Flow :	5 L/min
Injection Volume :	0.3	T.Flow :	0.2 ml/min
Sample Name :	KRAS WT Peptide (VVVGAGGVGK)	B.conc :	50%H2O/50%MeOH
Cat No. :	RP30235		
Lot No. :	E1975124A		
Theoretical MW:	842.00		
Observed MW :	841.8		