

Sample Name :OVA Peptide (323-339)

Cat No. :RP10610

Lot No. :E1993226F

Time Processed :12:33:38

Month-Day-Year Processed :06/08/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20510453

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.900	8029	1413	0.141
2	8.300	8242	1593	0.144
3	8.433	90699	13087	1.589
4	8.592	20813	4060	0.365
5	8.808	5425484	860149	95.065
6	9.008	105074	13509	1.841
7	9.242	6846	1644	0.120
8	9.458	41970	4679	0.735
Total		5707156	900134	100.000