

Sample Name : Influenza A NP(366-374) Strain A/PR/8/35
 Cat No. : RP20267
 Lot No. : E2004624G

Time Processed :15:43:52
 Month-Day-Year Processed :07/16/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
 Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	3.867	21991	4053	1.107
2	4.892	21476	1942	1.081
3	5.133	8751	1590	0.440
4	5.317	1928566	269737	97.054
5	5.633	6325	1002	0.318
Total		1987108	278323	100.000