

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E2038924J

Time Processed :17:55:37
Month-Day-Year Processed :10/26/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.583	56240	6021	0.939
2	11.108	5895364	581820	98.443
3	11.733	22100	1393	0.369
4	12.275	14876	1779	0.248
Total		5988580	591013	100.000