


Mass Spectrum



Sample Information

Month-Day Processed :	10/25/24	Equipment	:ZJ22010150
Time Processed :	16:05:07	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	T.Flow	:0.2 ml/min
Cat No. : RP11110		B.conc	:50%H2O/50%MeOH
Lot No. : E2038924J			
Theoretical MW : 2197.60			
Observed MW : 2197.2			