


Mass Spectrum



Sample Information

Month-Day Processed :	08/19/24	Equipment	:ZJ22010150
Time Processed :	13:39:47	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human		
Cat No. :	RP11110	T.Flow	:0.2 ml/min
Lot No. :	E2071024G	B.conc	:50%H2O/50%MeOH
Theoretical MW :	2197.60		
Observed MW :	2197.6		