

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E2143925B
Time Processed :6:56:13
Month-Day-Year Processed :02/18/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.917	23181	2404	0.223
2	20.150	37634	2848	0.362
3	20.508	35865	3620	0.345
4	20.775	29247	3216	0.281
5	21.058	60485	4788	0.581
6	21.458	10128589	970966	97.315
7	21.883	93016	9032	0.894
Total		10408018	996874	100.000