


# Mass Spectrum



## Sample Information

|                          |                    |                 |                 |
|--------------------------|--------------------|-----------------|-----------------|
| Month-Day Processed :    | 02/18/25           | Equipment       | :ZJ22010150     |
| Time Processed :         | 8:58:06            | Interface Bias  | : +4.5 kV       |
| Injection Volume :       | 0.2                | Drying Gas Flow | :5 L/min        |
| Sample Name :            | Melittin,honey bee | T.Flow          | :0.2 ml/min     |
| Cat No. : RP20415        |                    | B.conc          | :50%H2O/50%MeOH |
| Lot No. : E2143925B      |                    |                 |                 |
| Theoretical MW : 2847.50 |                    |                 |                 |
| Observed MW : 2847.6     |                    |                 |                 |