


Mass Spectrum



Sample Information

Month-Day Processed : 02/13/23
Time Processed : 15:38:43
Injection Volume : 0.4
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E2309123A

Theoretical MW:1621.91

Observed MW : 1621.6

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment
Interface Bias
Drying Gas Flow
T.Flow
B.conc

: ZJ21010035
: +4.5 kV
:5 L/min
:0.2 ml/min
:50%H₂O/50%MeOH