

Sample Name :KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E2312925A

Time Processed :11:03:09

Month-Day-Year Processed :01/20/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010140

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.808	4508	1080	0.141
2	7.942	66451	10148	2.077
3	8.258	4989	726	0.156
4	8.833	3287	677	0.103
5	9.142	3040767	356969	95.059
6	9.425	11525	2355	0.360
7	10.308	43915	5771	1.373
8	10.425	7541	1478	0.236
9	12.975	4928	734	0.154
10	14.817	10905	1378	0.341
Total		3198815	381318	100.000