


Mass Spectrum



Sample Information

Month-Day Processed :	01/20/25	Equipment	: ZJ21010035
Time Processed :	10:50:02	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	KRAS G12V Peptide (VVVGAVGVGK)	T.Flow	: 0.2 ml/min
Cat No. :	RP30229	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E2312925A		
Theoretical MW :	884.08		
Observed MW :	884.2		