

Sample Name :TAT peptide
Cat No. : RP20256
Lot No. : E2314523F
Time Processed :17:32:39
Month-Day-Year Processed :06/04/2023
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	0
5.00	Pumps	Pump A B.Conc	0
20.00	Pumps	Pump A B.Conc	30
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	0
35.00	Pumps	Pump A B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK11010011

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.760	9646	1150	0.093
2	19.998	38136	3755	0.369
3	20.175	13485	2305	0.131
4	20.292	17911	2981	0.173
5	20.490	9872398	744958	95.605
6	21.264	162663	7813	1.575
7	21.924	42190	2715	0.409
8	22.567	17032	1353	0.165
9	22.756	47429	5195	0.459
10	23.007	87786	9079	0.850
11	23.474	8140	1185	0.079
12	24.018	9438	1005	0.091
Total		10326252	783492	100.000