

Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat

Cat No. : RP11090

Lot No. : E2379426D

Time Processed : 1:26:20

Month-Day-Year Processed : 04/09/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010053

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.558	18160	1836	0.317
2	12.858	5676463	737147	98.947
3	13.417	16638	1610	0.290
4	13.658	25634	2306	0.447
Total		5736895	742899	100.000