


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	08/13/26	Interface	: ZI21010035
Time Processed :	13:27:46	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	TAT peptide	T.Flow	: 0.2 ml/min
Cat No. : RP20256		B. conc	: 50% H ₂ O/50% MeOH
Lot No. : E2416026H			
Theoretical MW :	1621.91		
Observed MW :	1622.0		