

Sample Name :OVA Peptide (323-339)

Cat No. :RP10610

Lot No. :E2439126F

Time Processed :14:14:28

Month-Day-Year Processed :06/26/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ24010222

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.108	4050	824	0.074
2	7.283	38029	7065	0.699
3	7.650	5181068	811993	95.297
4	7.817	186083	31588	3.423
5	8.217	14485	1816	0.266
6	8.308	9151	1634	0.168
7	8.417	3871	1000	0.071
Total		5436736	855920	100.000