

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E2440726F
Time Processed :19:49:58
Month-Day-Year Processed :06/25/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ25010054

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.567	11590	1710	0.191
2	7.717	74834	8732	1.234
3	8.092	5783381	848993	95.345
4	8.275	156157	29552	2.574
5	8.608	15874	2920	0.262
6	8.650	23917	3024	0.394
Total		6065753	894931	100.000