

Sample Name :KRAS G12C Peptide (VVVGACGVGK)

Cat No. :RP30236

Lot No. :E2564626F

Time Processed :19:41:58

Month-Day-Year Processed :06/20/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.758	23390	2921	0.392
2	9.125	112725	10383	1.889
3	9.883	26004	5996	0.436
4	10.050	5720671	651831	95.866
5	11.167	76904	12799	1.289
6	11.458	7677	1573	0.129
Total		5967371	685502	100.000