

Sample Name:OVA Peptide (323-339)

Cat No. : RP10610

Lot No. : E2565626F

Time Processed :11:32:31

Month-Day-Year Processed :06/17/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.877	42004	7171	0.182
2	8.237	531807	64930	2.304
3	8.481	379446	62257	1.644
4	8.740	21929137	2239367	95.000
5	9.162	164979	20195	0.715
6	9.361	22117	4067	0.096
7	9.591	13871	3358	0.060
Total		23083360	2401347	100.000