


Mass Spectrum



Sample Information

Month-Day Processed :	12/09/24	Equipment	:ZJ22010150
Time Processed :	20:16:27	Interface Bias	:+4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	β-Amyloid (25-35)	T.Flow	:0.2 ml/min
Cat No. :	RP10008	B.conc	:50%H2O/50%MeOH
Lot No. :	E2569324L		
Theoretical MW:	1060.27		
Observed MW :	1060.0		