

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E2577224L

Time Processed :14:05:50

Month-Day-Year Processed :12/11/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010508

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.333	12502	1741	0.605
2	17.783	21695	2893	1.050
3	18.167	25976	3184	1.257
4	18.958	1975944	282324	95.588
5	19.292	26460	2702	1.280
6	19.467	4574	953	0.221
Total		2067151	293795	100.000