


Mass Spectrum



Sample Information

Month-Day Processed : 12/09/24
 Time Processed : 20:16:08
 Injection Volume : 0.2
 Sample Name : Exendin-4
 Cat No. : RP10874
 Lot No. : E2581824L
 Theoretical MW: 4186.66
 Observed MW : 4186.8

Interface

Nebulizing Gas Flow :1.5L/min
 CDL Temp :250
 Block Temp :200

:ESI

Equipment : ZJ21010035
 Interface Bias : +4.5 kV
 Drying Gas Flow :5 L/min
 T.Flow :0.2 ml/min
 B.conc :50%H₂O/50%MeOH