

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human

Cat No. :RP11110

Lot No. :E2584525L

Time Processed :23:28:18

Month-Day-Year Processed :12/21/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010348

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	15.050	49043	6823	0.766
2	15.383	26755	3303	0.418
3	15.650	27301	5419	0.427
4	15.833	6096355	839144	95.242
5	16.000	201474	62806	3.148
Total		6400928	917495	100.000