


Mass Spectrum



Sample Information		Interface	:ESI	Equipment	:ZJ22010150
Month-Day Processed :	12/21/25	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Time Processed :	18:44:30	CDL Temp	:250	Drying Gas Flow	:5 L/min
Injection Volume :	0.3	Block Temp	:200	T.Flow	:0.2 ml/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human				
Cat No. :	RP11110	B.conc			:50%H2O/50%MeOH
Lot No. :	E2584525L				
Theoretical MW :	2197.60				
Observed MW :	2197.8				