


Mass Spectrum



Sample Information

Month-Day Processed : 12/10/24

Time Processed : 8:58:01

Injection Volume : 0.4

Sample Name : TAT peptide

Cat No. : RP20256

Lot No. : E2593824L

Theoretical MW : 1621.91

Observed MW : 1622.0

Interface :ESI

Nebulizing Gas Flow :1.5L/min

CDL Temp :250

Block Temp :200

Equipment : ZJ21010035

Interface Bias : +4.5 kV

Drying Gas Flow :5 L/min

T.Flow :0.2 ml/min

B.conc :50%H2O/50%MeOH