

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E2595224L

Time Processed :9:38:18

Month-Day-Year Processed :12/21/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0310

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.833	270220	10949	2.471
2	11.350	11183	2125	0.102
3	11.642	60122	10690	0.550
4	11.883	4560	840	0.042
5	12.225	10441050	1282912	95.459
6	12.458	123319	24202	1.127
7	14.142	27310	2688	0.250
Total		10937765	1334406	100.000