


Mass Spectrum



Sample Information		Equipment : ZJ21010035	
Month-Day Processed :	12/11/24	Interface Bias :	+4.5 kV
Time Processed :	9:25:57	Drying Gas Flow :	5 L/min
Injection Volume :	0.2	T.Flow :	0.2 ml/min
Sample Name :	RGD	B.conc :	50%H2O/50%MeOH
Cat No. :	RP20297		
Lot No. :	E2597224L		
Theoretical MW :	1025.06		
Observed MW :	1025.0		

Interface :	ESI
Nebulizing Gas Flow :	1.5L/min
CDL Temp :	250
Block Temp :	200