


Mass Spectrum



Sample Information

Month-Day Processed :	04/27/23	Equipment	:ZJ22010150
Time Processed :	7:40:25	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	:5 L/min
Sample Name :	β-Amyloid (25-35)	T.Flow	:0.2 ml/min
Cat No. :	RP10008	B.conc	:50%H2O/50%MeOH
Lot No. :	E2651523D		
Theoretical MW:1060.27			
Observed MW : 1060.2			