


Mass Spectrum



Sample Information

Month-Day Processed :	06/28/23	Equipment :	ZJ21010035
Time Processed :	9:39:24	Interface Bias :	+4.5 kV
Injection Volume :	0.2	Drying Gas Flow :	5 L/min
Sample Name :	RGD	T.Flow :	0.2 ml/min
Cat No. :	RP20297	B.conc	50%H2O/50%MeOH
Lot No. :	E2654623F		
Theoretical MW :	1025.06		
Observed MW :	1025.0		