


# Mass Spectrum



## Sample Information

Month-Day Processed :	01/14/25	Equipment	:ZJ22010150
Time Processed :	22:49:32	Interface Bias	:+4.5 kV
Injection Volume :	0.4	Drying Gas Flow	:5 L/min
Sample Name :	DYKDDDDDK Peptide	T.Flow	:0.2 ml/min
Cat No. :	RP10586	B.conc	:50%H2O/50%MeOH
Lot No. :	E2659725A		
Theoretical MW :	1012.98		
Observed MW :	1012.6		