


Mass Spectrum



Sample Information

Month-Day Processed :	11/08/24	Equipment	:ZJ22010150
Time Processed :	10:50:54	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	DYKDDDDK Peptide	T.Flow	:0.2 ml/min
Cat No. : RP10586		B.conc	:50%H2O/50%MeOH
Lot No. : E2680124K			
Theoretical MW: 1012.97			
Observed MW : 1012.6			