

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E2852824I

Time Processed :16:33:21

Month-Day-Year Processed :10/05/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
10.00	Pumps	B.Conc	35
35.00	Pumps	B.Conc	95
37.00	Pumps	B.Conc	95
37.01	Pumps	B.Conc	5
45.00	Pumps	B.Conc	5
45.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	15.865	8280	842	0.127
2	16.241	30969	2738	0.477
3	16.360	10704	1970	0.165
4	16.565	51523	5192	0.793
5	16.904	18866	1936	0.290
6	17.088	38080	7186	0.586
7	17.304	6221567	818691	95.769
8	17.563	98043	12033	1.509
9	18.004	766	-6	0.012
10	18.870	17665	2661	0.272
Total		6496463	853242	100.000