

Sample Name: Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E2859224I

Time Processed: 15:22:01

Month-Day-Year Processed :09/28/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment:ZJ19010014

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.195	44326	7139	1.993
2	8.032	18458	3251	0.830
3	9.841	2114954	335798	95.100
4	10.013	46180	6457	2.077
Total		2223918	352645	100.000