

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. :RP10769

Lot No. :E2952825J

Time Processed :15:38:40

Month-Day-Year Processed :10/26/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ23010112

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.692	73767	2089	0.931
2	16.775	32309	3514	0.408
3	16.975	115535	7692	1.458
4	17.317	70145	6552	0.885
5	17.617	7615981	1182201	96.087
6	17.900	18358	2555	0.232
Total		7926095	1204603	100.000