


Mass Spectrum



Sample Information

Month-Day Processed : 03/08/26
Time Processed : 17:26:15
Injection Volume : 0.3
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E2972026C
Theoretical MW: 1621.91
Observed MW: 1622.0

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment :ZJ22010150
Interface Bias :+4.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H2O/50%MeOH