

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E2979025C
Time Processed :20:16:56
Month-Day-Year Processed :04/07/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm


<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ17010509

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm				
Peak#	Ret. Time	Area	Height	Area%
1	12.667	9515	1085	0.381
2	13.158	2481234	338638	99.362
3	13.483	6413	1147	0.257
Total		2497162	340871	100.000