


# Mass Spectrum



## Sample Information

Month-Day Processed : 04/07/25  
Time Processed : 21:43:02  
Injection Volume : 0.3  
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human  
Cat No. : RP11110  
Lot No. : E2979025C  
Theoretical MW : 2197.60  
Observed MW : 2197.5

Equipment : ZJ21010035  
Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH