


Mass Spectrum



Sample Information

Month-Day Processed :	05/18/22	Equipment	: ZJ21010035
Time Processed :	16:31:19	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	TAT peptide	T.Flow	: 0.2 ml/min
Cat No. :	RP20256	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E3004122E		
Theoretical MW :	1621.91		
Observed MW :	1621.6		