


Mass Spectrum



Sample Information

Month-Day Processed : 02/19/24
 Time Processed : 8:38:22 PM
 Injection Volume : 0.3
 Sample Name : Exendin-4
 Cat No. : RP10874
 Lot No. : E3039824A
 Theoretical MW: 4186.66
 Observed MW : 4186.4

Interface

Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

:ESI

Equipment : ZJ22010150
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH