


Mass Spectrum



Sample Information

Month-Day Processed : 12/02/23

Time Processed : 8:53:46 AM

Injection Volume : 0.4

Sample Name : KRAS WT Peptide (VVVGAGGVGK)

Cat No. : RP30235

Lot No. : E3040523L

Theoretical MW : 842.00

Observed MW : 841.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ22010150
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH