

Sample Name:[Gly22]- β -Amyloid (1-42), Arctic Mutation
Cat No. : RP30247
Lot No. : E3208124L
Time Processed :13:36:55
Month-Day-Year Processed :12/08/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :UniPSN 30-300 30um,4.6x250mm

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.883	104445	2162	0.939
2	15.925	109975	6972	0.989
3	16.958	10692927	289353	96.154
4	18.242	150464	4624	1.353
5	18.992	62815	1996	0.565
Total		11120626	305107	100.000