


Mass Spectrum



Sample Information

Month-Day Processed : 07/26/26

Time Processed : 11:00:01

Injection Volume : 0.4

Sample Name : C-Peptide (3-33), human

Cat No. : RP11116

Lot No. : E3257726G

Theoretical MW : 3020.26

Observed MW : 3019.6

Interface

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment

Equipment : ZJ22010150

Interface Bias : -3.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50%H₂O/50%MeOH