


# Mass Spectrum



Sample Information  
Month-Day Processed : 12/19/25  
Time Processed : 2:59:38 PM  
Injection Volume : 0.4  
Sample Name : SARS-CoV-2 Spike-derived Peptide (YLQPRTFLL)  
Cat No. : RP30230  
Lot No. : E3467525L  
Theoretical MW: 1150.38  
Observed MW : 1150.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : ZJ22010150  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH