

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. :RP17313
 Lot No. :E3493321L
 Time Processed :12:08:25
 Month-Day-Year Processed :12/31/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.325	2041	342	0.006
2	11.125	2452	260	0.007
3	11.461	65168	8549	0.186
4	11.700	13625	1204	0.039
5	12.204	241664	22013	0.691
6	12.404	33722644	1640641	96.469
7	12.829	804456	84540	2.301
8	13.367	16301	1604	0.047
9	13.629	7079	942	0.020
10	13.850	4543	747	0.013
11	14.312	77120	5809	0.221
Total		34957093	1766650	100.000