

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E3507924J

Time Processed :19:58:26

Month-Day-Year Processed :11/03/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
10.00	Pumps	B.Conc	35
35.00	Pumps	B.Conc	95
37.00	Pumps	B.Conc	95
37.01	Pumps	B.Conc	5
45.00	Pumps	B.Conc	5
45.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ23010113

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.544	13133	376	0.285
2	13.568	11771	325	0.255
3	13.975	4235	774	0.092
4	15.219	2801	307	0.061
5	16.249	2831	378	0.061
6	17.117	4554763	634778	98.683
7	19.150	4210	742	0.091
8	22.873	2113	248	0.046
9	25.305	8563	914	0.186
10	33.450	11126	615	0.241
Total		4615544	639457	100.000