


Mass Spectrum



Sample Information

Month-Day Processed :	10/31/24	Interface	:ESI	Equipment	:ZJ22010150
Time Processed :	22:35:32	Nebulizing Gas Flow	:1.5L/min	Interface Bias	:+4.5 kV
Injection Volume :	0.3	CDL Temp	:250	Drying Gas Flow	:5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	Block Temp	:200	T.Flow	:0.2 ml/min
Cat No. :	RP10769	B.conc		B.conc	:50%H2O/50%MeOH
Lot No. :	E3507924J				
Theoretical MW :	3766.20				
Observed MW :	3766.2				