

Sample Name: Tau peptide (275-305) (Repeat2 domain)
 Cat No. : RP30887
 Lot No. : E3572826C
 Time Processed : 0:07:19
 Month-Day-Year Processed : 03/28/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm
 <<LC Time Program>>


Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm
 Equipment: ZJ19010142

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.113	13162	4272	0.720
2	10.291	1775058	282172	97.059
3	10.521	12341	2808	0.675
4	11.273	15795	3323	0.864
5	18.405	12483	2274	0.683
Total		1828840	294848	100.000