


Mass Spectrum



Sample Information

Month-Day Processed : 10/27/24

Time Processed : 0:29:51

Injection Volume : 0.4

Sample Name : DYKDDDDDK Peptide

Cat No. : RP10586

Lot No. : E3642424J

Theoretical MW : 1012.97

Observed MW : 1012.8

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment : ZJ21010035

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50% H<sub>2</sub>O / 50% MeOH