

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. : E3732326B
Time Processed :3:46:27
Month-Day-Year Processed :02/23/2026
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.758	13026	2192	0.294
2	8.892	49549	8281	1.118
3	9.083	4309139	715833	97.192
4	9.350	14016	4631	0.316
5	9.525	8246	1642	0.186
6	12.850	30949	4334	0.698
7	13.325	8714	1665	0.197
Total		4433640	738578	100.000